

ProBelt – Perfectly hygienic CT belt and fixation systems

ProBelt is the first completely hook-and-loop-free belt system and allows a perfectly hygienic fixation of patients on the CT. The ProBelt system is compatible with the CT tables of all common equipment manufacturers using different rail adapters.



No hook-and-loop



Perfectly hygienic

- Smooth PU-foil without hook-and-loop fastener
- Safe to disinfect even with 90% alcohol



High Durability


- Durable and hard-wearing materials with high-quality workmanship
- Strong immobilization even with uncooperative patients



Quick and easy to use

- Rail adapter to move ProBelt to the right position
- Simple and intuitive closing mechanism





ProBelt in Practice

Dr. Sascha Pietruschka, Chief Medical Officer Institute of Radiology and Neuroradiology, DBK-Neubrandenburg

"[...] It is out of the question to avoid a contamination of the conventional 'fabric-based belts/straps' provided by the CT-equipment manufacturers. At the same time, a quick and proper disinfection bringing down the contamination of these straps below the maximum tolerable level, is also impossible. [...]"

We were able to prove that by using commercially available disinfectants, a hygienically perfect and immediate reuse of the 'ProBelt' is possible and that the "ProBelt" system solves the hygiene problem 100 %.

"[...] According to my radiographers, the handling and fixation strength is at least as good as with conventional belts."

Yes, I would like to learn more about the ProBelt fixation system

Please contact me for:

Product demonstration on-site

Additional information

A non-binding offer

First name _____

Function _____

Name _____

Phone _____

Company _____

E-mail _____

Please send this form or an e-mail to: info@pearl-technology.ch